


RUO

GASTROENTERITIS RT-qPCR MX-24 PANEL



Gastroenteritis infections may cause public health emergencies. It is essential to investigate the source and responsible pathogen causing gastroenteritis infections. **Bio-Speedy® Gastroenteritis RT-qPCR MX-24 Panel** helps to identify 24 pathogens from research samples.



10 SAMPLES / **65 MINUTES**

VIRUS

Sapovirus (GI/GII/GIV/GV)
Adenovirus
Astrovirus
Norovirus (GI/GII)
Rotavirus (A)

PARASITE

Giardia lamblia
Entamoeba histolytica
Cryptosporidium spp.
Cyclospora cayetanensis

GRAM POSITIVE BACTERIA

Clostridium difficile toxin A
Clostridium difficile toxin B
Clostridium difficile Binary toxin A/B

GRAM NEGATIVE BACTERIA

Plesiomonas shigelloides
Salmonella spp.
Campylobacter spp.
Vibrio parahaemolyticus
Vibrio cholerae
Shigella/Enteroinvasive *E. coli* (EIEC)
Enterohaggagative *E. coli* (EAEC)
Shiga-toxin producing *E. coli* (STEC)
Enteropathogenic *E. coli* (EPEC)
Enterotoxigenic *E. coli* (ETEC)
Yersinia enterocolitica

Compatible Instruments: Bio-Rad CFX96 Touch, Bio-Rad CFX Opus 96, Bio-Rad CFX384, Bio-Rad CFX Opus 384, Bioer LineGene 9600 Plus, Bioer QuantGene 9600 Real Time PCR Systems, Quant Studio 5, Quant Studio 6, Quant Studio 7

PRODUCT	PRODUCT SIZE	CATALOG NO
Bio-Speedy® Gastroenteritis RT-qPCR MX-24S Panel	100 Rxns	BS-GE-MX24S-100R
Bio-Speedy® Gastroenteritis RT-qPCR MX-24T Panel	100 Rxns	BS-GE-MX24T-100R

For Research Use Only (RUO). Not for use in diagnostic procedures. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of disease.



US: 23232 Peralta Dr. STE 220 Laguna Hills CA 92653
Main: Huzur Mah. Metin Oktay Cad. NuroI Life Sitesi D Blok No:3/31 Sariyer-Istanbul-TURKEY
Phone: (949) 685-7834 - www.bioeksen-us.com - e-mail: info@axisbio.org