


STI MINI RT-qPCR KIT



STDs are infections transmitted among people through sexual contact, and the detection of the pathogen is crucial in order to prevent transmission, provide treatment and reduce the risk of negative health effects. *Mycoplasma genitalium*, *Neisseria gonorrhoeae*, *Chlamydia trachomatis* and *Trichomonas vaginalis* are the most prevalent bacteria causing STI. **Bio-Speedy® STI Mini RT-qPCR Kit** targets these 4 pathogens and provide comprehensive insights for diagnosis.



94 SAMPLES / **65 MINUTES**

Compatible Instruments: Bio-Rad CFX96 Touch, Bio-Rad CFX Opus 96, Bioer LineGene 9600 Plus, Bioer QuantGene 9600 Real Time PCR Systems

GRAM POSITIVE BACTERIA

Neisseria gonorrhoeae
Trichomonas vaginalis

GRAM NEGATIVE BACTERIA

Mycoplasma genitalium
Chlamydia trachomatis

PRODUCT	PRODUCT SIZE	CATALOG NO
Bio-Speedy® STI Mini RT-qPCR Kit	100 Rxns	BS-SM-100R

For Research Use Only (RUO). Not for use in diagnostic procedures. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of disease.