


RUO

# SEXUALLY TRANSMITTED INFECTIONS RT-qPCR PANEL



STDs are infections transmitted among people by sexual contact. Pathogen studies are crucial in order to understand the mechanism and finding out the way to diagnosis and treatment. **Bio-Speedy® Sexually Transmitted Infections RT-qPCR Panel** covers 13 STD related pathogens.



22 SAMPLES / **65** MINUTES

**Compatible Instruments:** Bio-Rad CFX96 Touch, Bio-Rad CFX Opus 96, Bio-Rad CFX384, Bio-Rad CFX Opus 384, Bioer LineGene 9600 Plus, Bioer QuantGene 9600 Real Time PCR Systems, Quant Studio 5, Quant Studio 6, Quant Studio 7

## VIRUS

HSV-1

HSV-2

## GRAM POSITIVE BACTERIA

*Gardnerella vaginalis*

*Neisseria gonorrhoeae*

*Treponema pallidum*

*Trichomonas vaginalis*

*Ureaplasma parvum/urealyticum*

*Streptococcus agalactiae*

## GRAM NEGATIVE BACTERIA

*Mycoplasma genitalium*

*Mycoplasma hominis*

*Chlamydia trachomatis*

*Haemophilus ducreyi*

PRODUCT	PRODUCT SIZE	CATALOG NO
Bio-Speedy® Sexually Transmitted Infections RT-qPCR Panel	100 Rxns	BS-STI-S-100R
Bio-Speedy® Sexually Transmitted Infections RT-qPCR Panel	100 Rxns	BS-STI-T-100R

For Research Use Only (RUO). Not for use in diagnostic procedures. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of disease.