



SIGMOIDA ANALYSIS SOFTWARE

- Automated result interpretation
- Reporting in 2 seconds
- The **offline** software ensures data safety
- Automated report transfer to the Laboratory Information Systems.
- Human error elimination
- Ready-to-use for multiple kits



Mic qPCR Cycler Plate Setup Interface

Test Kit : Bio-Speedy® Kit
Samples : 40
Positive Control : 4
Negative Control : 4

	FAM	HEX	ROX	Cy5	1	2	3	4	5	6	7	8	9	10	11	12
A	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC
B	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC
C	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC
D	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC
E	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC
F	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC
G	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC
H	S	S	S	S	S	S	S	S	S	S	S	S	S	S	NC	PC

PC : Positive Control NC : Negative Control Fluorophores : FAM / HEX / ROX / Cy5

Bio-Rad Plate Setup Interface

Report

Test Kit : Bio-Speedy® Kit
File : CVD FLU DATA 1.pcrd
Date : 03 Feb 2023 09:10:22
User : Admin
Version : V 2.1.3

Show in plate view

All (171) Positive (37) Negative (101) Issues (33) Controls (39)

Select All

	Well	Sample	Target	Cq	Result	Warnings
<input checked="" type="checkbox"/>	A1	CVD HASTA 1	SARS-CoV-2	22.32	Positive	
<input checked="" type="checkbox"/>	A1	CVD HASTA 1	Influenza B	N/A	Negative	
<input checked="" type="checkbox"/>	A1	CVD HASTA 1	Influenza A	N/A	Negative	
<input checked="" type="checkbox"/>	B1	CVD HASTA 1	SARS-CoV-2	24.04	Positive	
<input checked="" type="checkbox"/>	B1	CVD HASTA 1	Influenza B	N/A	Negative	
<input checked="" type="checkbox"/>	B1	CVD HASTA 1	Influenza A	N/A	Negative	
<input checked="" type="checkbox"/>	D1	INF A HASTA 1	SARS-CoV-2	N/A	Negative	
<input checked="" type="checkbox"/>	D1	INF A HASTA 1	Influenza B	N/A	Negative	
<input checked="" type="checkbox"/>	D1	INF A HASTA 1	Influenza A	N/A	Negative	
<input checked="" type="checkbox"/>	E1	INF B HASTA 1	SARS-CoV-2	N/A	Negative	
<input checked="" type="checkbox"/>	E1	INF B HASTA 1	Influenza B	N/A	Negative	
<input checked="" type="checkbox"/>	E1	INF B HASTA 1	Influenza A	N/A	Negative	
<input checked="" type="checkbox"/>	A2	CVD HASTA 2	SARS-CoV-2	25.2	Positive	
<input checked="" type="checkbox"/>	A2	CVD HASTA 2	Influenza B	N/A	Negative	
<input checked="" type="checkbox"/>	A2	CVD HASTA 2	Influenza A	N/A	Negative	
<input checked="" type="checkbox"/>	B2	CVD HASTA 2	SARS-CoV-2	24.81	Positive	
<input checked="" type="checkbox"/>	B2	CVD HASTA 2	Influenza B	N/A	Negative	
<input checked="" type="checkbox"/>	B2	CVD HASTA 2	Influenza A	N/A	Negative	

Result Reporting Table